

MONIDOR

Precise intravenous therapy for every patient

Monidrop[®] W

INFUSION MONITOR

Enables more precise intravenous therapy and helps nurses to improve patient safety

High quality

Small and accurate device

More precise IV therapy

In order to improve patient safety

User friendly

Easy to set up, train and service

Cost-saving

Suitable for standard IV administration sets that are up to 80 % lower price



IV ScreenTM

Better patient management and quicker response times

Overview

IV Screen will show the status of all the infusions on the hospital wards

Efficient

Saves nurses time as they are able to monitor all infusions remotely

Less alarms

Notifications are given remotely which is especially beneficial at night-time



# WSOEDN							Sim Ward 🚽 📃
7 Room 1 1	15%	6 Room 1 2 NaCl	99% III I	8 Room 2 1 NaCl	97% 880	5 Room 2 2 NaCl	85% III
Ot 58min		3t 56min		3t Omin	A Ot 13min / 1t Omin	1t Omin	
n 240 ml/h		n 95 ml/h		🕋 80 ml/h	₩ 100 ml/h	n 50 ml/h	
218 ml		루 395 ml		루 223 ml	🔺 41 ml / 100 ml	루 50 ml	
2 Room 3 2	925 BB	3 Room 4 1	97% 💷	9 Room 4 2	995 000	4 Room 4 3	205.000
1t 28min		1t 16min		Ot 58min		1t ómin	
n 100 ml/h		n 195 ml/h		🕋 47 ml/h		🕋 210 ml/h	
루 146 ml		루 284 ml		📮 43 ml		👎 218 ml	
10 Room 5 1	92N 💷	1 Room 5 2	65% BID				
🕒 1t 46min		🕒 1t 8min					
🥋 125 ml/h		n 220 ml/h					
寻 211 ml		📮 280 ml					
			tronilitaatie 3, 90590 Oulu, F Tattoruoja B	enseenstelle jappens eely RCAD-1-103 32 195 5021 (elegeneerke on researchers			
			tronilitaatie 3, 90590 Oulu, F Tattoruoja B	9NLAND +358 50 295 9063 info@monitor.co Inkitamolota			
			tronilitaatie 3, 90590 Oulu, F Tattoruoja B	9NLAND +358 50 295 9063 info@monitor.co Inkitamolota			
			tronilitaatie 3, 90590 Oulu, F Tattoruoja B	9NLAND +358 50 295 9063 info@monitor.co Inkitamolota		_	
-			tronilitaatie 3, 90590 Oulu, F Tattoruoja B	9NLAND +358 50 295 9063 info@monitor.co Inkitamolota			
			tronilitaatie 3, 90590 Oulu, F Tattoruoja B	9NLAND +358 50 295 9063 info@monitor.co Inkitamolota	-		



IV Screen[™]

Overview screen shows

- Status of all infusions
- IV Screen shows remotely:
 - 1. Notification status (traffic light)
 - 2. Flow rate (ml/h)
 - 3. Volume (ml)
 - 4. Duration of treatment (hh:min)
 - 5. Battery status
 - 6. Target settings (if set)

MONIDOR

For availability please contact:

Monidor Oy Elektroniikkatie 3 90590 Oulu Finland

+358 102959063 info@monidor.com www.monidor.com

